



Genome polyprotein (2039-2279)

OVERVIEW

PRODUCT REFERENCE	PX-P4343
ORIGIN SPECIES	
FRAGMENT	Glu2039-Gln2279
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	55.36kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P04517
ALIAS/SYNONYMS	Mature chain, Nuclear inclusion protein A

