

N-acyl homoserine lactone hydrolase (SAMN04487767_102296)

OVERVIEW

PRODUCT REFERENCE	PX-P4359
ORIGIN SPECIES	
FRAGMENT	Met1-Ile250
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	29.20KDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	A0A1D3P3H1
ALIAS/SYNONYMS	N-acyl homoserine lactonase/N-acyl homoserine lactone hydrolase