



Antiviral protein DAP-30

OVERVIEW

PRODUCT REFERENCE	PX-P4372
ORIGIN SPECIES	
FRAGMENT	Ala24-Lys277
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	29.84kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P24476 (99,2%)
ALIAS/SYNONYMS	Dianthin 30,Ribosome-inactivating protein,rRNA N-glycosidase

