



DNA-binding protein H-NS (hns)

OVERVIEW

PRODUCT REFERENCE	PX-P4380
ORIGIN SPECIES	
FRAGMENT	Met12-Gln156
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	17.63kDa
FORM	
BUFFER	20 mM Tris HCl, pH 8.0, 2 mM MgCl ₂ , and 20 mM NaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P0ACF9
ALIAS/SYNONYMS	DNA-binding protein H-NS, Histone-like protein HLP-II, Protein B1, Protein H1

