



DsRed protein (DsRed)

OVERVIEW

PRODUCT REFERENCE	PX-P4388
ORIGIN SPECIES	
FRAGMENT	Leu3-Arg230
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	31.03KDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q2WG74
ALIAS/SYNONYMS	DsRed/DsRed protein

