



Peripheral plasma membrane protein CASK(CASK) (1-926)

OVERVIEW

PRODUCT REFERENCE	PX-P4393
ORIGIN SPECIES	
FRAGMENT	Met1-Try926
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	106KDa
FORM	
BUFFER	50 mM Tris-HCl pH 7.5, 150 mM NaCl, 10% Glycerol
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O14936
ALIAS/SYNONYMS	CASK/Peripheral plasma membrane protein CASK

