



Heme sensor protein HssS(hssS)

OVERVIEW

PRODUCT REFERENCE	PX-P4404
ORIGIN SPECIES	
FRAGMENT	Ile188-Leu457
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	40.73KDa
FORM	
BUFFER	TBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	A8Z553
ALIAS/SYNONYMS	hssS/Heme sensor protein HssS

