



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Listeriolysin O(hly)

OVERVIEW

PRODUCT REFERENCE	PX-P4407
ORIGIN SPECIES	
FRAGMENT	Asp26-Glu529
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	56.94KDa
FORM	
BUFFER	50 mM Tris-HCl pH 8.0, 150 mM NaCl, 3mM Cysteine, 0.3mM Cystine
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P13128
ALIAS/SYNONYMS	hly/Listeriolysin O



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.