



Proteogenix

Antibodies and Proteins

Proteogenix SAS
Espace Européen de l'Entreprise
19 rue de La Haye
67300 Schiltigheim, France
Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.Proteogenix.science

ORF256(orf256)

OVERVIEW

CATALOG NO	PX-P4409
ORIGIN SPECIES	
FRAGMENT	Cys54-Gly256
MOLECULAR WEIGHT	24.73KDa
DELIVERED WITH TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Procaryotic expression
HOST	Escherichia coli (E.coli)
PURITY	90% by SDS-PAGE
UNIPROT ID	Q7GI75
FORM	Liquid
BUFFER	50 mM Tris-HCl pH 7.5, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

More info

SYNONYMS

orf256/ORF256

UNIPROT LINK

<http://www.uniprot.org/uni-prot/Q7GI75>

SEQUENCE

MGSHHHHHHSGCNLFGASFIITLPPEIQDPQALAHLAGLNFYLSLYEQDPGWVTFIQNELNHNTPLEDIPGRLKFLMEEKLSSMRQDVIQEF
VALYQRIGPYLPIEPYLVEALRSYLDHIHATDSFTVLQASYQDLRENEGGSVFFRDAVSHNRDLLEAESSARRCLEVEQRIRWEEIPKSKAS
LERAHEHALDLFKSEDLRRELEKKRAG