



Chromosomal replication initiator protein DnaA

OVERVIEW

PRODUCT REFERENCE	PX-P4506
ORIGIN SPECIES	
FRAGMENT	Mer1-472Ser
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	54.43kDa
FORM	
BUFFER	50 mM Tris-HCl pH 8.0, 150 mM NaCl/PBS pH 7.5, 4 M urea
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	D7HH42
ALIAS/SYNONYMS	/

