



Phosphatidylinositol 4-kinase alpha 1(PI4KA1) (1-344)

OVERVIEW

PRODUCT REFERENCE	PX-P4509
ORIGIN SPECIES	
FRAGMENT	Met1-Gln344
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	38.83kDa
FORM	
BUFFER	PBS pH 7.5, 4 M urea
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9SXA1
ALIAS/SYNONYMS	Phosphatidylinositol 4-kinase alpha 1, PI4-kinase alpha 1, PtdIns-4-kinase alpha 1, Phosphatidylinositol 4-OH kinase alpha1, AtPI4Kalpha1, PI-4Kalpha1

