



Ribosomal RNA small subunit methyltransferase H(rsmH)

OVERVIEW

PRODUCT REFERENCE	PX-P4517
ORIGIN SPECIES	
FRAGMENT	full length
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	35.96kDa
FORM	
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	B4TXH0
ALIAS/SYNONYMS	16S rRNA m(4)C1402 methyltransferase, rRNA (cytosine-N(4)-)-methyltransferase RsmH

