



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Alanine racemase, catabolic(dadX)

OVERVIEW

PRODUCT REFERENCE	PX-P4523
ORIGIN SPECIES	
FRAGMENT	full length
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	41.97kDa
FORM	
BUFFER	PBS, pH 7.5, 4M urea
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P06191
ALIAS/SYNONYMS	



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.