



Otogelin-like protein(Otogl)

OVERVIEW

| | |
|---------------------------|-----------------|
| PRODUCT REFERENCE | PX-P4536 |
| ORIGIN SPECIES | |
| FRAGMENT | Arg1507-Lys1727 |
| TAG | |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Prokaryotic expression |
| HOST | Escherichia coli (E.coli) |
| MOLECULAR WEIGHT | 25.90kDa |
| FORM | |
| BUFFER | PBS, pH7.5, 4M urea |
| STABILITY & STORAGE | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| UNIPROT | F7A4A7 |
| ALIAS/SYNONYMS | |

