

Lmo1842 protein(lmo1842)

OVERVIEW

PRODUCT REFERENCE	PX-P4542
ORIGIN SPECIES	
FRAGMENT	full length
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	32.43kDa
FORM	
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8Y658
ALIAS/SYNONYMS	

