



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Troponin C, slow skeletal and cardiac muscles(TNNC1)&Troponin I, cardiac muscle(TNNI3)

OVERVIEW

PRODUCT REFERENCE	PX-P4547
ORIGIN SPECIES	
FRAGMENT	Met1~Glu161 & Thr31~Arg111
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	28.85kDa
FORM	
BUFFER	50 mM Tris-HCl pH 7.5, 150 mM NaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P63316&P19429
ALIAS/SYNONYMS	TNNC&TNNC1,TN-C & Cardiac troponin I



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.