



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Rat Resistin-like alpha(Retnla)

OVERVIEW

PRODUCT REFERENCE	PX-P4554
ORIGIN SPECIES	Rattus norvegicus (Rat)
FRAGMENT	Met1~Ser111
TAG	C-terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	35-40kDa
FORM	
BUFFER	Tris-glycine
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q99P85
ALIAS/SYNONYMS	Fizz1, Cysteine-rich secreted protein FIZZ1, RELMalpha



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.