



Parathyroid hormone(PTH)

OVERVIEW

PRODUCT REFERENCE	PX-P4568
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Gln115
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	13.83kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P01270
ALIAS/SYNONYMS	Parathormone;Parathyrin

