



Tumor-associated calcium signal transducer 2(TACSTD2)

OVERVIEW

PRODUCT REFERENCE	PX-P4569
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Thr274
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	35-45kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P09758
ALIAS/SYNONYMS	GA733-1, M1S1, TROP2, Cell surface glycoprotein Trop-2 ; Membrane component chromosome 1 surface marker 1 ; Pancreatic carcinoma marker protein GA733-1

