



Interleukin-12 (Il12a & Il12b heterodimer) (Met1-Ala215)&(Met1-Ser335)

OVERVIEW

PRODUCT REFERENCE	PX-P4570
ORIGIN SPECIES	House Mouse
FRAGMENT	Il12a(Met1-Ala215)& Il12b(Met1-Ser335)
TAG	C-terminal His Tag & C-terminal Twin-Strep tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	Il12a(25kDa) & Il12b (45-60kDa)
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P43431&P43432
ALIAS/SYNONYMS	IL-12 subunit p35 & IL-12 subunit p40

