



Thioredoxin 1(trxA)(Shigella flexneri)

OVERVIEW

PRODUCT REFERENCE	PX-P4573
ORIGIN SPECIES	Shigella flexneri
FRAGMENT	Met1-Ala109
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	12.87kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P0AA30
ALIAS/SYNONYMS	Trx-1

