



Vascular endothelial Growth Factor proteins A(VEGFA)

OVERVIEW

PRODUCT REFERENCE	PX-P4574
ORIGIN SPECIES	
FRAGMENT	Met1-Arg191
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	20.13kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P15692-4
ALIAS/SYNONYMS	VEGF, VPF

