



Histone H3.1(HIST1H3A)

OVERVIEW

PRODUCT REFERENCE	PX-P4575
ORIGIN SPECIES	General
FRAGMENT	Met1~Ala136
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	16.47kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P68431/Q6LED0/Q6LBF0/P68433
ALIAS/SYNONYMS	H3FJ

