



Interleukin-12 (IL12AB heterodimer) (Met1-Ser219) &(Met1-Ser328)

OVERVIEW

PRODUCT REFERENCE	PX-P4576
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	II12A(Met1-Ser219) & II12B(Met1-Ser328)
TAG	C-terminal Strep tag & C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	40kDa & 43kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P29459&P29460
ALIAS/SYNONYMS	IL12AB

