

Homeobox protein engrailed-1 (EN1)

OVERVIEW

PRODUCT REFERENCE	PX-P4579
ORIGIN SPECIES	
FRAGMENT	Met1-Glu392
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	42.14kDa
FORM	
BUFFER	PO4(pH:7.4) 20 mM with about 0.5M NaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q05925
ALIAS/SYNONYMS	Homeobox protein en-1, Hu-En-1

