



Programmed cell death 1 ligand 1(Cd274)

OVERVIEW

PRODUCT REFERENCE	PX-P4580
ORIGIN SPECIES	House Mouse
FRAGMENT	Met1-His239
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	40-45kDa(26.92 kDa)
FORM	Lyophilized
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9EP73
ALIAS/SYNONYMS	B7h1, Pdcd1l1, Pdcd1lg1, Pdl1,PD-L1,PDCD1 ligand 1,Programmed death ligand 1,B7 homolog 1,B7-H1,CD274

