



Thy-1 membrane glycoprotein(THY1)

OVERVIEW

PRODUCT REFERENCE	PX-P4582
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Cys130
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	22-30kDa(14.55 kDa)
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P04216
ALIAS/SYNONYMS	CDw90,Thy-1 antigen, CD90

