

JAK2 protein - Tyrosine-protein kinase JAK2(JAK2)

OVERVIEW

PRODUCT REFERENCE	PX-P4590
ORIGIN SPECIES	
FRAGMENT	Asn508~Ala800
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	35kDa and 30kDa(35.93kDa)
FORM	
BUFFER	PBS pH 7.5, 0.02% SKL
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O60674
ALIAS/SYNONYMS	Janus kinase 2 , JAK-2

