



Peroxisome proliferator-activated receptor alpha(PPARA)

OVERVIEW

PRODUCT REFERENCE	PX-P4597
ORIGIN SPECIES	
FRAGMENT	Ala201-Tyr468
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	31.6kDa
FORM	
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q07869
ALIAS/SYNONYMS	Nuclear receptor subfamily 1 group C member 1, NR1C1, PPAR

