



N-acyl homoserine lactonase(aiiA)

OVERVIEW

PRODUCT REFERENCE	PX-P4602
ORIGIN SPECIES	
FRAGMENT	Met1-Ile250
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	29.3kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9L8R8
ALIAS/SYNONYMS	AHL-lactonase,Acyl-homoserine lactonase

