



Tubuliny-Tyr carboxypeptidase 2 (VASH2)

OVERVIEW

PRODUCT REFERENCE	PX-P4617
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Thr2-Ile354
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q86V25
ALIAS/SYNONYMS	Vasohibin-2, Vasohibin-like protein, VASHL

