



Sialidase-3(NEU3)

OVERVIEW

PRODUCT REFERENCE	PX-P4627
ORIGIN SPECIES	
FRAGMENT	Met1-Asn461
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9UQ49-2
ALIAS/SYNONYMS	Ganglioside sialidase,Membrane sialidase,N-acetyl-alpha-neuraminidase 3

