

Vitamin K-dependent protein C (PROC)

OVERVIEW

PRODUCT REFERENCE	PX-P4631
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Pro461
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P04070
ALIAS/SYNONYMS	Anticoagulant protein C, Autoprothrombin IIA, Blood coagulation factor XIV