



Leucine-rich repeats and immunoglobulin-like domains protein 2(LRIG2)

OVERVIEW

PRODUCT REFERENCE	PX-P4632
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Leu42-Thr805
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O94898
ALIAS/SYNONYMS	LIG-2,KIAA0806, LIG2

