



Interleukin-1 receptor-associated kinase 4 (Irak4)

OVERVIEW

PRODUCT REFERENCE	PX-P4633
ORIGIN SPECIES	House Mouse
FRAGMENT	Met1-Ala459
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8R4K2
ALIAS/SYNONYMS	IRAK-4

