



Interleukin-2 (IL2)

OVERVIEW

PRODUCT REFERENCE	PX-P4637
ORIGIN SPECIES	Oryctolagus cuniculus (Rabbit)
FRAGMENT	Ala21-Ser153
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O77620
ALIAS/SYNONYMS	IL-2, T-cell Growth Factor proteins, TCGF

