



Fibrinogen-like protein 1(FGL1)

OVERVIEW

PRODUCT REFERENCE	PX-P4638
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ile312
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	35.32kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q08830
ALIAS/SYNONYMS	HP-041,Hepassocin,HPS,Hepatocyte-derived fibrinogen-related protein 1,HFREP-1,Liver fibrinogen-related protein 1,LFIRE-1,HFREP1

