



T-cell surface glycoprotein CD8 alpha chain (CD8A)

OVERVIEW

PRODUCT REFERENCE	PX-P4639
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Asp182
TAG	C-terminal Twin-Strep Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P01732
ALIAS/SYNONYMS	T-lymphocyte differentiation antigen T8/Leu-2,CD8a,MAL

