

Immunoglobulin G-binding protein G(spg) His Tag

OVERVIEW

PRODUCT REFERENCE	PX-P4645
ORIGIN SPECIES	Streptococcus sp. group G
FRAGMENT	Lys314-Glu497
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	21.24kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P19909
ALIAS/SYNONYMS	IgG-binding protein G