



## Zinc transporter 8(SLC30A8)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4648
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Met275-Asp369 + PSTPPGSSGGG + Leu267-Asp369
<b>TAG</b>	C-terminal 10*His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Yeast
<b>MOLECULAR WEIGHT</b>	
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q8IWU4
<b>ALIAS/SYNONYMS</b>	ZnT-8, Solute carrier family 30 member 8,ZNT8

