



## Progranulin(Grn)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4649
<b>ORIGIN SPECIES</b>	
<b>FRAGMENT</b>	Met1-Leu589
<b>TAG</b>	C-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Insect
<b>MOLECULAR WEIGHT</b>	64.4kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P28798
<b>ALIAS/SYNONYMS</b>	Acrogranin, Epithelin/granulin precursor, Glycoprotein of 88 Kda, GP88, Glycoprotein 88, PC cell-derived Growth Factor proteins, PCDGF, Proepithelin, PEPI

