



Dual specificity protein phosphatase CDC14A(Cdc14a)

OVERVIEW

| | |
|---------------------------|--------------------|
| PRODUCT REFERENCE | PX-P4660 |
| ORIGIN SPECIES | House Mouse |
| FRAGMENT | Met1-Tyr603 |
| TAG | C-terminal His Tag |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Eukaryotic expression |
| HOST | Mammalian cells |
| MOLECULAR WEIGHT | |
| FORM | |
| BUFFER | PBS, pH7.5 |
| STABILITY & STORAGE | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| UNIPROT | Q6GQT0 |
| ALIAS/SYNONYMS | CDC14 cell division cycle 14 homolog A |

