



Claudin-3(CLDN3)

OVERVIEW

PRODUCT REFERENCE	PX-P4669
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Val220
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O15551
ALIAS/SYNONYMS	Clostridium perfringens enterotoxin receptor 2,CPE-R 2,CPE-receptor 2,Rat ventral prostate.1 protein homolog,hRVP1,C7orf1, CPETR2

