



Interleukin-27 subunit beta(EBI3)&Interleukin-12 subunit alpha(IL12A)

OVERVIEW

PRODUCT REFERENCE	PX-P4678
ORIGIN SPECIES	
FRAGMENT	Met1-Lys229(Q14213) & Arg23-Ser219(P29459)
TAG	C-terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	74.4kDa
FORM	Lyophilized
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q14213&P29459
ALIAS/SYNONYMS	IL-27 subunit beta,IL-27B,Epstein-Barr virus-induced gene 3 protein,EBV-induced gene 3 protein,IL27B & IL-12A,Cytotoxic lymphocyte maturation factor 35 kDa subunit,CLMF p35,IL-12 subunit p35,NK cell stimulatory factor chain 1,NKSF1,NKSF1

