



Thymic stromal lymphopoietin(TSLP)

OVERVIEW

PRODUCT REFERENCE	PX-P4679
ORIGIN SPECIES	
FRAGMENT	Tyr29-Gln159
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	18.2kDa
FORM	Lyophilized
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q969D9
ALIAS/SYNONYMS	Thymic stromal lymphopoietin

