



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Hyaluronan mediated motility receptor (HMMR)

OVERVIEW

PRODUCT REFERENCE	PX-P4687
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Lys651
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O75330-2
ALIAS/SYNONYMS	Intracellular hyaluronic acid-binding protein, Receptor for hyaluronan-mediated motility, CD168, IHABP, RHAMM



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.