



Runt-related transcription factor 1 (RUNX1)

OVERVIEW

PRODUCT REFERENCE	PX-P4688
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Tyr480
TAG	N-terminal His Tag & sumo Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	
MOLECULAR WEIGHT	64.8kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q01196-8
ALIAS/SYNONYMS	Acute myeloid leukemia 1 protein,Core-binding factor subunit alpha-2,CBF-alpha-2,Oncogene AML-1,Polyomavirus enhancer-binding protein 2 alpha B subunit,PEA2-alpha B,PEBP2-alpha B,SL3-3 enhancer factor 1 alpha B subunit,SL3/AKV core-binding factor alpha B

