



Interleukin-12 subunit alpha (IL12a) & Interleukin-12 subunit beta(IL12b)

OVERVIEW

| | |
|---------------------------|--|
| PRODUCT REFERENCE | PX-P4690 |
| ORIGIN SPECIES | House Mouse |
| FRAGMENT | Met1-Ala215 & Met1-Ser335 |
| TAG | C-terminal His Tag & C-terminal Twin-Strep Tag |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Eukaryotic expression |
| HOST | Mammalian cells |
| MOLECULAR WEIGHT | 40kDa & 43kDa |
| FORM | |
| BUFFER | PBS, pH7.5 |
| STABILITY & STORAGE | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| UNIPROT | P43431&P43432 |
| ALIAS/SYNONYMS | IL-12A,IL-12 subunit p35 & IL-12B,IL-12 subunit p40 |

