

Alpha-lactalbumin(LALBA)

OVERVIEW

PRODUCT REFERENCE	PX-P4700
ORIGIN SPECIES	
FRAGMENT	Lys20-Leu142
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	15.1kDa
FORM	
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P00709
ALIAS/SYNONYMS	LYZL7, Lactose synthase B protein, Lysozyme-like protein 7

