



RAC-beta serine, threonine-protein kinase(AKT2)

OVERVIEW

PRODUCT REFERENCE	PX-P4706
ORIGIN SPECIES	
FRAGMENT	Met1-Lys146
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	16.96 kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P31751
ALIAS/SYNONYMS	Protein kinase Akt-2, Protein kinase B beta, PKB beta, RAC-PK-beta

